

Cheat Sheet

For Clinicians



Calm practitioner equals calm patient

Step 1

Remove NeedleCalm™ device from the fridge.

Step 2

Prepare the procedure as per the instructions for use.

Step 3

Place, Press and Push.

Place NeedleCalm™ on the patient above the position of the intended needle procedure location and press on the NeedleCalm™ device as needle is pushed.

Needle placement and suggested hand position for:

(1) Intra-muscular injection, (2) Subcutaneous injection and (3) Intravenous Cannulation.

